




Copyright © Taiwan Surgical Corporation

Reusable Clip Applier

Reusable Clip Applier - Multi-fire with
single-use cartridge and reusable
handle 10/20 Titanium Clips, ML/L
10mm.

Outstanding control and durability

www.twsc.com.tw



Yellow Bar
Remind surgeon
when using out clips




Ring
Block interstitial fluid
getting into handle
during replacement

Anti-glare Shaft
Minimize reflection and
instrument moving during firing



Reusable Handle
(#TRC-H0000)



Ergonomic Rotation knob
The angle is adjustable during
firing process



Single-use Cartridge, 20 Clips, ML 10mm (#TRC-TM121)



Single-use Cartridge, 20 Clips, L 10mm (#TRC-TL121)



Specifications:

Product No.	Description	Shaft Spec. (D*L)	Packaging
TRC-TM110	Titanium Clip Cartridge, 10 Clips, ML 10mm	ø10*309mm	25pcs / box, 4 boxes / carton
TRC-TM121	Titanium Clip Cartridge, 20 Clips, ML 10mm	ø10*309mm	25pcs / box, 4 boxes / carton
TRC-TL110	Titanium Clip Cartridge, 10 Clips, L 10mm	ø10*359mm	25pcs / box, 4 boxes / carton
TRC-TL121	Titanium Clip Cartridge, 20 Clips, L 10mm	ø10*359mm	25pcs / box, 4 boxes / carton
TRC-H0000	Reusable Clip Applier Handle		1 pcs / box, 24 boxes / carton